



Product catalogue POINT OF CARE



.....• pursuing the excellence



XLAB SOLUTIONS



XLAB SOLUTIONS

We are X-Lab Solutions

The story begins in 1992 when "a few good men" have founded a service organization for laboratory and blood gas instruments in Transilvania.

At first, X-Lab Solutions' activity focused on maintenance and repairs. Along the way, the development of turnkey medical projects and providing lab supplies, have become priorities.

Today, the main focus is on providing after-sales services for medical equipment and projects (installation, commissioning, training, application support, service, maintenance).

X-Lab Solutions achievements:

- OVER 110 equipped hospitals and the associated services delivered, all over Europe
- Over 120 private laboratories equipped with diagnostic systems and serviced
- Over 170 private healthcare institutions equipped with medical devices and accessories



XLAB SOLUTIONS


The pillars of X-Lab Solutions



POC
Point of Care



W'sHC
In Vitro Fertilisation
& Medical Devices



SP
Special Projects

Point of Care brings important benefits for patients and the medical team: it improves their experience through accurate results, reduced waiting time and increased availability.

Point of Care easily fulfills its mission of providing patient safety right at the bedside.



XLAB SOLUTIONS

Point Of Care products




Blood gas and
electrolyte analyzers



Cardiac markers
analyzers, sepsis,
pregnancy, coagulation



Transcutaneous
monitors



Collection
systems



Single parameter
POC analyzers



Multi parameter
POC analyzers



Urine
analyzers



Nurse-patient
communication systems



Laryngoscopes

POINT OF CARE

Urine analyzers



.....● pursuing the excellence



XLAB SOLUTIONS



XLAB SOLUTIONS



Laura & Laura Smart Semi-Automated & Compact Urine Strip Readers

Urine Analysis: Quick, Simple, Smart



Laura
Semi-Automated
Urine Strip Reader



Laura Smart
Compact
Urine Strip Reader



BLOOD GAS AND
ELECTROLYTE ANALYZERS



CARDIAC MARKERS ANALYZERS,
SEPSIS, PREGNANCY, COAGULATION



TRANSCUTANEOUS
MONITORS



COLLECTION
SYSTEMS



SINGLE PARAMETER
POC ANALYZERS



MULTI-PARAMETER
POC ANALYZERS



URINE
ANALYZERS

- Laura & Laura Smart



NURSE-PATIENT
COMMUNICATION SYSTEMS



LARYNGOSCOPES



About •

Urine strip readers LAURA & LAURA Smart

- Efficient, powerful and compact solutions for routine chemical urinalysis
- Semi-automated technology
- Facilitate the samples processing and contribute to correct diagnosis of the patients

Specifications:

LAURA

Power source: AC 90 – 230 V/50 – 60 Hz
Power consumption: 45 W max./6 W in standby mode
Temperature: 20 – 25 °C
Relative humidity: 20 – 80 %
Weight: 5 kg including the adaptor
Dimensions: 430 x 290 x 170 mm

LAURA Smart

Power source: Normal mode: 9 V DC switch adapter, AC 90 – 230 V/50 – 60 Hz
Portable mode: 6 x 1.5 V (type AA) batteries
Power consumption: 20 W max./1 W standby mode
Temperature: 20 – 25 °C
Relative humidity: 20 – 80 %
Weight: 0.8 kg without adapter, battery and paper roll
Dimensions: 230 x 127 x 110 mm

CAT. NR.	PRODUCT	QTY.
50001727	LAURA®	1 pce
10008297	DekaPHAN® LAURA	100 strips
10008298	HeptaPHAN® LAURA	100 strips
10010266	Urinorm	2 x 3 bottles
50003491	Control strips for LAURA®	5 pcs

CAT. NR.	PRODUCT	QTY.
50003508	LAURA® Smart	1 pce
10008297	DekaPHAN® LAURA	100 strips
10008298	HeptaPHAN® LAURA	100 strips
10010239	PentaPHAN® LAURA	100 strips
10010238	DiaPHAN® LAURA	100 strips
10010262	MicroalbuPHAN® LAURA	50 strips
10010266	Urinorm	2 x 3 bottles
50003513	Control strips for LAURA® Smart	3 pcs

 BLOOD GAS AND ELECTROLYTE ANALYZERS

 CARDIAC MARKERS ANALYZERS, SEPSIS, PREGNANCY, COAGULATION

 TRANSCUTANEOUS MONITORS

 COLLECTION SYSTEMS

 SINGLE PARAMETER POC ANALYZERS

 MULTI-PARAMETER POC ANALYZERS

 URINE ANALYZERS

● Laura & Laura Smart

 NURSE-PATIENT COMMUNICATION SYSTEMS

 LARYNGOSCOPES



Advantages •

Urine strip readers LAURA & LAURA Smart


 BLOOD GAS AND ELECTROLYTE ANALYZERS

 CARDIAC MARKERS ANALYZERS, SEPSIS, PREGNANCY, COAGULATION

 TRANSCUTANEOUS MONITORS

 COLLECTION SYSTEMS

 SINGLE PARAMETER POC ANALYZERS

 MULTI-PARAMETER POC ANALYZERS

 URINE ANALYZERS

• Laura & Laura Smart

 NURSE-PATIENT COMMUNICATION SYSTEMS

 LARYNGOSCOPES

Control, calibration and maintenance

Simple and fast automatic self-check and calibration at each start-up

Flexibility

Compatible with different types of urine test strips

Two and three sample identification modes

LAURA Smart: Sequential, Sample ID
LAURA: Sequential, Sample ID and Worklist

Compatibility with Erba Lachema and 3rd party liquid controls



Total availability in less than 24 h anywhere in Romania and Moldova.

Contact us for **technical support, service for your equipment or consumables.**



We are offering an **unique and customized range of services:**

SALES SERVICES

Our sales representatives deeply understand your needs and guide in taking your best decisions.

TRAINING SERVICES

Through our training services you will know exactly how to operate your equipment according to the manufacturer's instructions.

MAINTENANCE & SUPPORT

Our engineers are highly skilled in servicing the equipment we provide.

Cluj-Napoca Head Office

Calea Manastur 70, Ap. 30, Cluj-Napoca, Romania
Tel.: +40 264 42 76 77 • E-mail: office@xlab.ro

Bucharest

Sos. Odai 207A, Bucharest 1, Romania
Tel.: +40 21 332 76 31 • E-mail: office@xlab.ro

Chisinau

Negresteni 9, Chisinau, Moldova
Tel.: +40 729 038 316 • E-mail: office@xlab.ro

www.xlab.ro