



Micropipettes & Strippers W`SHC



XLAB SOLUTIONS

.....• pursuing the excellence



XLAB SOLUTIONS

We are X-Lab Solutions

The story begins in 1992 when "a few good men" have founded a service organization for laboratory and blood gas instruments in Transilvania.

At first, X-Lab Solutions' activity focused on maintenance and repairs. Along the way, the development of turnkey medical projects and providing lab supplies, have become priorities.

Today, the main focus is on providing after-sales services for medical equipment and projects (installation, commissioning, training, application support, service, maintenance).

X-Lab Solutions achievements:

- OVER 110 equipped hospitals and the associated services delivered, all over Europe
- Over 120 private laboratories equipped with diagnostic systems and serviced
- Over 170 private healthcare institutions equipped with medical devices and accessories



XLAB SOLUTIONS


The pillars of X-Lab Solutions



POC
Point of care



W`sHC
In Vitro Fertilisation
& Medical Devices



SP
Special Projects

W`sHC - delivering excellence at every stage of the journey to a healthy and happy life!

IVF offers perfectly controlled conditions and state-of-the-art equipment to make every woman's dream of becoming a mother come true.

In Vitro Fertilisation

Micropipettes & Strippers



About •—————

ICSI pipettes

Designed to meet the most specific preferences and optimise results.

TPC Sharp-Blunt Spike™ Pipette - LISR

Designed for precise fluid control and suction.

Easy penetration and gentle progression through the oolemma, minimizing the risk of tearing.



ORIGIO® SlimLine Pipette – MIC-SLM-30

Short taper for added stability when injecting.

Smaller inner diameter to minimize the impact on the oocyte.



TPC Pipette



ORIGIO® Pipette



-  ICSI Pipettes
-  Holding Pipettes
-  Biopsy Pipettes
-  ORIGIO® Micropipettes
-  TPC™ Micropipettes
-  The Stripper® & EZ-Grip®

About ●

Holding pipettes


Designed for stable embryo handling during biopsy and micromanipulation.

-  ICSI Pipettes
-  Holding Pipettes
-  Biopsy Pipettes
-  ORIGIO[®] Micropipettes
-  TPC[™] Micropipettes
-  The Stripper[®] & EZ-Grip[®]




Small

Small: Inner diameter
15-20 μm / Outer Diameter
65-95 μm



Medium

Medium: 15-20 μm /
Outer diameter
95-120 μm




Large

Large: Inner diameter
25-30 μm /Outer
Diameter 120-150 μm



Extra Large










XL: Inner Diameter
25-30 μm /Outer
diameter 150-180 μm




About •—————

Biopsy pipettes

Different sizes to cover most procedures, from polar body to blastomere or trophoctoderm biopsies.

Inner Diameter	Design options
 <p>Polar Body</p> <p>14 μm</p>	 <p>FP - Flat & polished</p>  <p>BP - beveled & polished</p>  <p>BS - beveled & spiked</p>
 <p>Blastomere</p> <p>30, 35 or 40 μm</p>	 <p>FP - Flat & polished</p>  <p>BP - beveled & polished</p>
 <p>Trophoctoderm</p> <p>20 or 25 μm</p>	 <p>FP - Flat & polished</p>



-  ICSI Pipettes
-  Holding Pipettes
-  Biopsy Pipettes
-  ORIGIO® Micropipettes
-  TPC™ Micropipettes
-  The Stripper® & EZ-Grip®



XLAB SOLUTIONS



ICSI Pipettes



Holding Pipettes



Biopsy Pipettes




ORIGIO[®] Micropipettes



TPC[™] Micropipettes



The Stripper[®] & EZ-Grip[®]



CooperSurgical[®]

Fertility Solutions

ORIGIO[®] Micropipettes





XLAB SOLUTIONS



ICSI Pipettes



Holding Pipettes



Biopsy Pipettes




ORIGIO® Micropipettes



TPC™ Micropipettes



The Stripper® & EZ-Grip®



CooperSurgical®

Fertility Solutions

TPC™ Micropipettes


Quality-assured design by embryologists for embryologists





XLAB SOLUTIONS

-  ICSI Pipettes
-  Holding Pipettes
-  Biopsy Pipettes
-  ORIGIO® Micropipettes
-  TPC™ Micropipettes
-  The Stripper® & EZ-Grip®



CooperSurgical®
Fertility Solutions

The Stripper® & EZ-Grip®

Tools for denudation and handling





XLAB SOLUTIONS

- ICSI Pipettes
- Holding Pipettes
- Biopsy Pipettes
- ORIGIO® Micropipettes
- TPC™ Micropipettes
- The Stripper® & EZ-Grip®

About • _____

The Stripper® & EZ-Grip®

Adjustable

Aspiration volumes can easily be changed when needed

Compatible

Fits most plastic denudation and handling tips (0.9 mm diameter)

Proven

Tried and tested

Durable

Ergonomic design requiring minimal maintenance

Fast identification

EZ-Grip features a pipette size tip indicator dial

Varied

Pipettors to accommodate all sizes of pipettes

- Close-tolerance plastic tips are flexible, with smooth edges and consistent geometries
- Made from flexible medical-grade plastic, ensuring no scratches on Petri dishes
- Also available in a multi-pack vial for cost-effective options
- Designed for optimal performance and economy





XLAB SOLUTIONS

-  ICSI Pipettes
-  Holding Pipettes
-  Biopsy Pipettes
-  ORIGIO® Micropipettes
-  TPC™ Micropipettes
-  The Stripper® & EZ-Grip®

Bulb pipetting range

EZ-Squeeze™

The EZ-Squeeze is a medical grade plastic tip fitted with a flexible bulb for a reliable seal and simpler set-up. EZ-Squeeze tips are available in all common sizes.

EZ-Squeeze™ Handle

The autoclavable EZ-Squeeze handle is designed to offer ease of handling and a comfortable working position.



EZ-Strip®

Developed for strict lab sterility, minimizing cross-contamination risk. The long version is ideal for the handling of oocytes in tubes.



Total availability in less than 24 h anywhere in Romania and Moldova.

Contact us for technical support, service for your equipment or consumables.



We are offering an unique and customized range of services:

SALES SERVICES

Our sales representatives deeply understand your needs and guide in taking your best decisions.

TRAINING SERVICES

Through our training services you will know exactly how to operate your equipment according to the manufacturer's instructions.

MAINTENANCE & SUPPORT

Our engineers are highly skilled in servicing the equipment we provide.

Cluj-Napoca Head Office

Calea Manastur 70, Ap 30, Cluj-Napoca, Romania
Tel.: +40 264 42 76 77 · E-mail: office@xlab.ro

Bucharest

Sos. Odai 207A Bucharest 1, Romania
Tel.: +40 21 332 76 31 · E-mail: office@xlab.ro

Chisinau

Negresteni 9, Chisinau, Moldova
Tel.: +40 729038 316 · E-mail: office@xlab.ro

www.xlab.ro